

Soluzione Totale Genius 8080

Numero dell'articolo: XX-XXXX



Descrizione

Il Genius 1053 e il Genius 3070 sono costruiti attorno alla comprovata tecnologia Peaks, fornendo una soluzione a gas totale per il triplo quadrupolo Shimadzu LCMS-8080.



Generator Details

- **Tipo di gas:** Azoto
- **Compressore integrato:** No

Caratteristiche principali

- **Flessibile** - Disponibile con o senza compressore d'aria interno
- **Purezza** - Zero Air con meno di 0,1 ppm di contenuto di idrocarburi (come metano)
- **Durevole** - Poche parti mobili con requisiti minimi di manutenzione
- **Soluzione Single Source** - Un generatore può soddisfare i requisiti di più applicazioni
- **Installazione semplice** - Generatore progettato come plug & sistema di gioco
- **Economico** - Più conveniente rispetto a qualsiasi altro metodo di fornitura di gas
- **Conveniente** - Gas su richiesta, nessun rischio per la salute, non c'è bisogno di preoccuparsi di rimanere senza gas

Specifiche tecniche

Soluzione Totale Genius 8080	
Tipo di gas	Azoto
Tempo di avvio:	60 mins
Corrente:	0
Livello di rumore:	0

[PEAK Protected]TM

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: www.peakscientific.com/ordering/
- For service plans visit: www.peakscientific.com/service/service-plans/

Peak Scientific UK
Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8200

Peak Scientific North America
Tel: +1-800-767-6532
Fax: +1-978-608-9503

Peak Scientific China
Tel: +86 21 5079 1190
Fax: +86 21 5079 1191
For a full list of our worldwide office locations, please visit:
Web: www.peakscientific.com
Email: discover@peakscientific.com

Peak Scientific's Quality Management System conforms to: ISO:9001

