

# ZA070A Zero Air

Numero dell'articolo: UK: 30-1610  
(230v) | EU: 30-1620 (230v) | US: 30-  
1630 (110v)




## Descrizione

Generando aria zero per il gas del rivelatore nelle applicazioni GC, il generatore ZA070A rimuove tutti gli idrocarburi da una fonte di aria compressa prima che entrino nella strumentazione, assicurando i risultati più ordinati dal rivelatore.



## Generator Details

- **Tipo di gas:** Aria
- **Compressore integrato:** No



## Caratteristiche principali

- Produce un flusso costante di aria di grado zero con un contenuto di idrocarburi inferiore a 0,1 ppm
- Esegue 24 ore su 24, 7 giorni su 7, se necessario
- Un generatore può fornire più strumenti
- Poche parti in movimento significano poca manutenzione e garantiscono una lunga durata del generatore
- Richiede una configurazione minima
- Garanzia on-site completa di 12 mesi

## Specifiche tecniche

	ZA070A
Tipo di gas	Aria
Flusso massimo di gas:	7L/min
Contenuto di idrocarburi:	0.1ppm
Purezza massima:	0.1ppm THC
Tempo di avvio:	40 Minutes
Voltaggio:	230 ± 10% VAC
Frequenza:	50 HZ
Corrente:	4.9 Amps
Temperatura operativa massima:	30°C / 86°F
Particelle:	0.01µm
Dimensioni (AxLxP) in mm	880 x 430 x 410 mm
Dimensioni (AxLxP) in pollici	34,6 x 16,9 x 16,1 inches
Peso del generatore	64kg / 140,8lbs

# [PEAK Protected]<sup>TM</sup>

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: [www.peakscientific.com/ordering/](http://www.peakscientific.com/ordering/)
- For service plans visit: [www.peakscientific.com/service/service-plans/](http://www.peakscientific.com/service/service-plans/)

**Peak Scientific UK**  
Tel: +44 (0)141 812 8100  
Fax: +44 (0)141 812 8200

**Peak Scientific North America**  
Tel: +1-800-767-6532  
Fax: +1-978-608-9503

**Peak Scientific China**  
Tel: +86 21 5079 1190  
Fax: +86 21 5079 1191  
For a full list of our worldwide office locations, please visit:  
Web: [www.peakscientific.com](http://www.peakscientific.com)  
Email: [discover@peakscientific.com](mailto:discover@peakscientific.com)

Peak Scientific's Quality Management System conforms to: ISO:9001

